


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Bis(trimethylstannyl)acetylene

Product Number: 483613
CAS Number: 2117-50-2
MDL: MFCD01073798
Formula: C₈H₁₈Sn₂
Formula Weight: 351.65 g/mol



TEST

Specification

Appearance (Color)	White to Off-White
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Gravimetric Analysis	60.8 - 74.3 %
% Tin, After Ignition to Oxide	

Specification: PRD.0.ZQ5.10000055479

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.