

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Titanium(II) hydride - -325 mesh, 98%

**Product Number:** 209279  
**CAS Number:** 7704-98-5  
**MDL:** MFCD00151301  
**Formula:** H<sub>2</sub>Ti  
**Formula Weight:** 49.88 g/mol



TEST	Specification
Appearance (Color)	Grey to Very Dark Grey
Appearance (Form)	Powder
X-Ray Diffraction	Conforms to Structure
Gravimetric Analysis % Ti (After Ignition To Oxide)	93.5 - 98.3 %
Size -325 mesh	Conforms

Specification: PRD.1.ZQ5.10000017641

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.