

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Poly(tetrahydrofuran) – average $M_n \sim 2,000$, contains 200–350 ppm BHT as stabilizer

Product Number:

345326


CAS Number:

25190-06-1

MDL:

MFCD00148879

Formula:

 $H(OCH_2CH_2CH_2CH_2)_nOH$ **TEST****Specification**

Appearance (Color)

White

Appearance (Form)

Solid or Semi-solid

Infrared spectrum

Conforms to Structure

Molecular Weight

1900 - 2100

Color Test

 ≤ 40 APHA

Hydroxyl Value

53 - 59

Stabilizer

Confirmed

200-350ppm BHT

Specification: PRD.1.ZQ5.10000019798

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.