


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Melamine - 99%

Product Number: M2659
CAS Number: 108-78-1
MDL: MFCD00006055
Formula: C₃H₆N₆
Formula Weight: 126.12 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Purity (Titration by HCL04)	98.5 - 101.5 %
Solubility (Turbidity) 25mg/ml, H ₂ O	Clear to Very Slightly Hazy
Solubility (Color)	Colorless

Specification: PRD.1.ZQ5.10000030899

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.