

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Sodium hexametaphosphate - crystalline, +200 mesh, 96%

Product Number:

305553

CAS Number:

68915-31-1



MDL:

MFCD00131858

TEST	Specification
------	---------------

Appearance (Color)

White

Appearance (Form)

Crystals

ICP Major Analysis

Confirmed

Confirms Na and P Components

Purity (Titration by NaOH)

29.0 - 31.7 %

% P

Particle Size Distribution

+ 40 Mesh

Particle Size Distribution

+ 50 Mesh

Particle Size Distribution

+ 60 Mesh

Particle Size Distribution

+ 80 Mesh

Particle Size Distribution

+ 200 Mesh

Specification: PRD.0.ZQ5.10000039990

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.