

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium bromide - ReagentPlus®, ≥99%

Product Number: 220345
CAS Number: 7647-15-6
MDL: MFCD00003475
Formula: BrNa
Formula Weight: 102.89 g/mol

NaBr

| TEST | Specification |
|--|-----------------------|
| Appearance (Color) | White |
| Appearance (Form) | Powder or Crystals |
| X-Ray Diffraction | Conforms to Structure |
| Titration by AgNO ₃ | 99.0 - 101.0 % |
| Titration by AgNO ₃ % Br | 76.8 - 78.4 % |

Specification: PRD.0.ZQ5.10000016135

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.