

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

3-Allyloxy-2-hydroxy-1-propanesulfonic acid sodium salt solution - 40 wt. % in H₂O

Product Number:

409421

MDL:


MFCD00080396

Formula:

C₆H₁₁NaO₅S

Formula Weight:

218.2 g/mol

**TEST****Specification**

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Proton NMR spectrum

Conforms to Structure

Color Test

≤ 100 APHA

pH

6.0 - 9.0

10%; Water

Active Material

39.0 - 42.0 %

Sodium Bisulfite

≤ 0.1 %

Specification: PRD.1.ZQ5.10000009457.000

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.