

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Silver nanospheres – 200 nm avg. part. size, 0.02 mg/mL (in 2mM aqueous sodium citrate), citrate functionalized

Product Number: 796166

Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color)	Light Yellow to Yellow and Light Brown to Brown
Appearance (Form)	Liquid
Wavelength	
Report Result	
ICP Major Analysis	Confirmed
Confirms Silver Component	
Diameter	
Hydrodynamic (Report Result)	
Surface Area	
Report Result (TEM)	
Concentration	
Mass (Ag)(Report Result)	
Particles/ml	Report
Diameter	Confirmed
200 nm	
Miscellaneous Supplier Data	
Zeta Potential (mV)(Report Result)	
pH	
Report Result	
Silver (Ag)	Confirmed
99.99%	
Miscellaneous Supplier Data	Confirmed
Particle Surface: Sodium Citrate	
Recommended Retest Period	-----
12 months	

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Silver nanospheres – 200 nm avg. part. size, 0.02 mg/mL (in 2mM aqueous sodium citrate), citrate functionalized

Product Number: 796166

Storage Temperature: 2 - 8 °C

TEST	Specification
-------------	----------------------

Specification: PRD.2.ZQ5.10000098834

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.