

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

PAMAM dendrimer, ethylenediamine core, generation 1.0 solution - 20 wt. % in methanol

Product Number:

412384



MDL:

MFCD00192448

Formula:

C62H128N26O12

Formula Weight:

1,429.85 g/mol

Storage Temperature:

2 - 8 °C

TEST**Specification**

Appearance (Color)

Colorless to Light Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Residue on Evaporation

17.0 - 23.0 %

Specification: PRD.1.ZQ5.10000015400

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.