


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,1'-Bi-2-naphthol - 99%

Product Number: 104655
CAS Number: 602-09-5
MDL: MFCD00004068
Formula: C₂₀H₁₄O₂
Formula Weight: 286.32 g/mol



TEST

Specification

Appearance (Color)	White to Light Yellow and Faint Beige to Light Beige and Faint Brown to Light Brown
Appearance (Form)	Powder or Crystals or Chunk(s)
Infrared spectrum	Conforms to Structure
Purity (HPLC)	≥ 98.5 %
Purity (TLC)	≥ 99 %
Solubility (Turbidity)	Clear
0.5g/10ml, Dioxane	
Solubility (Color)	Colorless to Yellow

Specification: PRD.0.ZQ5.10000065594

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.