

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Ethyl 2-formyl-1-cyclopropanecarboxylate, predominantly trans - 96%

Product Number:**157279**

CAS Number:

20417-61-2

MDL:

MFCD00001285

Formula:


C₇H₁₀O₃

Formula Weight:

142.15 g/mol

Storage Temperature:

2 - 8 °C

**TEST****Specification**

Appearance (Color)

Colorless to Light Yellow

Appearance (Form)

Liquid

Proton NMR spectra

Conforms to Structure

Purity (GC)

≥ 95.5 %

Specification: PRD.1.ZQ5.10000005432

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.