

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Potassium - chunks (in mineral oil), 98% trace metals basis

Product Number:

244864

CAS Number:

7440-09-7

MDL:

MFCD00133776

Formula:

K

Formula Weight:

39.1 g/mol

**K**

TEST	Specification
------	---------------

Appearance (Color)

Conforms to Requirements

Grey to Dark Grey to Grey-Black

Appearance (Form)

Conforms

Sticks, Rods, or Chunks

Purity

Meets Requirements

98% Based on Trace Metals Analysis

Trace Metal Analysis

≤ 25000.0 ppm

Specification: PRD.2.ZQ5.10000005145

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.