

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Amino-1-propanethiol hydrochloride - technical grade

Product Number: **739294**
CAS Number: 7211-54-3



Formula: C₃H₉NS · HCl
Formula Weight: 127.64 g/mol

TEST	Specification
Appearance (Color)	White to Yellow
Appearance (Form)	Powder or Crystals
¹ H NMR Spectrum	Conforms to Structure
Titration with AgNO ₃ After Oxygen Combustion	80.0 - 120.0 %
Carbon	22.6 - 33.9 %
Nitrogen	8.8 - 13.2 %

Specification: PRD.3.ZQ5.10000015123

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.