

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Platinum(II) bromide – 98%

Product Number: 206121  
CAS Number: 13455-12-4  
MDL: MFCD00011180  
Formula: Br<sub>2</sub>Pt  
Formula Weight: 354.89 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Light Brown to Dark Brown
Appearance (Form)	Powder
ICP Major Analysis	Confirmed
Confirms Platinum Component	
Gravimetric Analysis	53.6 - 56.3 %
%Pt After Reduction With Hydrazine	

Specification: PRD.0.ZQ5.10000052036

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.