

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Hyaluronidase from bovine testes – Type I-S, lyophilized powder, 400–1000 units/mg solid

Product Number:

**H3506**

CAS Number:

37326-33-3

MDL:

MFCD00131351

Storage Temperature:

-20 °C

**TEST****Specification**

Appearance (Color)

Faint Beige to Beige to Brown

Appearance (Form)

Powder

units/mg solid

400 - 1000

One unit will cause a change in % transmittance at 600 nm of 0.330 per minute at pH 5.35 at 37 deg C in a 2.0 mL reaction mixture (45 minute assay).

Recommended Retest Period

2 years

Specification: PRD.3.ZQ5.10000015916

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.