

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Complement C3 from human serum - ≥ 3000 C3H50 units/mg (using C3 deficient serum)

Product Number:

C2910

MDL:

MFCD00130836

TEST	Specification
Protein Content	1.0 - 1.2 mg/ml
Purity by SDS Electrophoresis	$\geq 85\%$
units/mg protein	≥ 3000
Hemolytic titer C3H50	
HBSAG	Negative
HCV	Negative
HIV-1/HIV-2	Negative
Human Source Testing (Source)	Tested
Tested For Infectious Agents	Tested

Specification: PRD.1.ZQ5.10000006369

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.