

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Ytterbium(III) nitrate pentahydrate - 99.9% trace metals basis

**Product Number:****209147**Yb(NO<sub>3</sub>)<sub>3</sub> · 5H<sub>2</sub>O

CAS Number:

35725-34-9

MDL:

MFCD00149883

Formula:

N3O9Yb · 5H2O

Formula Weight:

449.13 g/mol

TEST	Specification
Appearance (Color)	White
Appearance (Form) Crystals or Crystalline Powder and/or Chunks	Conforms to Requirements
Complexometric EDTA % Yb	37.4 - 39.7 %
ICP Major Analysis Confirms Yb Component	Confirmed
Purity 99.9% Based On Rare Earth Analysis	Meets Requirements
Trace Rare Earth Analysis	≤ 1500.0 ppm

Specification: PRD.0.ZQ5.10000042679

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.