

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Carbonyldihydrido tris(triphenylphosphine)ruthenium(II)

Product Number: 335002
CAS Number: 25360-32-1
MDL: MFCD00015870
Formula: C₅₅H₄₇OP₃Rh
Formula Weight: 919.79 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White to Grey
Appearance (Form)	Conforms to Requirements
Powder and/or Lumps	
Infrared spectrum	Conforms to Structure
ICP Major Analysis	Confirmed
Confirms Ruthenium Component	
Carbon	69.4 - 74.5 %

Specification: PRD.1.ZQ5.10000028195

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.