


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3,5-Bis(trifluoromethyl)benzylamine - technical grade, 80%

Product Number: **263389**
CAS Number: 85068-29-7
MDL: MFCD00009909
Formula: C₉H₇F₆N
Formula Weight: 243.15 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Infrared spectrum
Purity (Titration by HCL04)
Purity (GC)

White to Yellow
Solid or Powder or Crystals or Chunks
Conforms to Structure
78.0 - 122.0 %
≥ 78.0 %

Specification: PRD.0.ZQ5.10000014903

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.