

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Monoclonal Anti-Human IgG2-FITC antibody produced in mouse - clone HP-6014, purified immunoglobulin, buffered aqueous solution

Product Number:

F4516

TEST	Specification
------	---------------

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Yellow
IgG Concentration (OD 280nm)	1.0 - 3.0 mg/ml
F/P Molar Ratio	3.0 - 8.0
Dot Blot-Titer	1:32-1:64
P-IMFA Particle Immunofl Assay	-----
F-PIMFA Agarose Beads Coating	Human IgG or Human IgG2
F-PIMFA Titer	1:16-1:32
Recommended Retest Period	-----
3 Years	

Specification: PRE.0.ZQ5.10000017057

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.