

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Diethylene glycol dimethyl ether - anhydrous, 99.5%

Product Number:

281662

CAS Number:

111-96-6

MDL:


MFCD00008503

Formula:

C₆H₁₄O₃

Formula Weight:

134.17 g/mol

**TEST****Specification**

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 99.45 %

Water (by Karl Fischer)

≤ 0.005 %

Residue on Evaporation

≤ 0.0005 %

Specification: PRD.0.ZQ5.10000020981

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.