


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium hydrogencitrate sesquihydrate - ReagentPlus[®], 99%

Product Number: 359084
CAS Number: 6132-05-4
MDL: MFCD00150445
Formula: C₆H₆Na₂O₇ · 1.5H₂O
Formula Weight: 263.11 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Purity (Titration by NaOH)	98.5 - 101.5 %
GC (area %)	≥ 98.5 %
Recommended Retest Period	-----
3 years	

Specification: PRD.0.ZQ5.10000078206

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.