

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Aluminum oxide – activated, basic, Brockmann I

Product Number:**199443**

CAS Number:

1344-28-1

MDL:

MFCD00003424

Formula:

Al₂O₃

Formula Weight:

101.96 g/mol



TEST	Specification
------	---------------

Appearance (Color)

White

Appearance (Form)

Powder

Vendor Information

Confirmed

Brockmann Number (Activity): I

pH Aqueous Suspension (Stirred): 9.5 +/- 0.5

Specific Surface Area: 205 m²/g +/- 10%Cl⁻: NILNa₂O: < 0.4%Fe₂O₃: < 0.03%SiO₂: < 0.03%

Specification: PRD.1.ZQ5.10000003722

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.