


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Bis(4-nitrophenyl) carbonate - $\geq 99\%$

Product Number: 161691
CAS Number: 5070-13-3
MDL: MFCD00007322
Formula: C₁₃H₈N₂O₇
Formula Weight: 304.21 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
White to Pale Yellow or Beige	
Appearance (Form)	Conforms to Requirements
Powder, Chips or Flakes and/or Chunks	
Infrared spectrum	Conforms to Structure
Carbon	50.8 - 51.8 %
Nitrogen	8.9 - 9.5 %

Specification: PRD.0.ZQ5.10000042248

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.