

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

4-Bromobenzylzinc bromide solution - 0.5 M in THF

Product Number:

499463

MDL:

MFCD01311435

Formula:


C₇H₆Br₂Zn

Formula Weight:

315.32 g/mol

Storage Temperature:

2 - 8 °C

**TEST****Specification**Purity Titration by EDTA
Complexometric

0.40 - 0.60 M

Specification: PRD.1.ZQ5.1000026099

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.