


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Methoxytyramine hydrochloride - ≥95.5%, crystalline

Product Number: M4251
CAS Number: 1477-68-5
MDL: MFCD00012902
Formula: C₉H₁₃NO₂ · HCl
Formula Weight: 203.67 g/mol
Storage Temperature: -20 °C



TEST

Specification

Appearance (Color)	White to Brown
Appearance (Form)	Conforms to Requirements
Crystalline Powder	
Infrared Spectrum	Conforms to Structure
Purity (Titration by HClO ₄)	≥ 95.5 %
Solubility (Turbidity)	Clear
50 mg/mL in Water	
Solubility (Color)	Colorless to Yellow to Brown
Recommended Retest Period	-----
5 Years	

Specification: PRD.0.ZQ5.10000089318

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.