

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Clostridium difficile Toxin A – lyophilized powder

Product Number:

C3977

MDL:

MFCD02098001

Storage Temperature:

2 - 8 °C

TEST**Specification**

Appearance (Form)

Lyophilized Powder

Purity by SDS Electrophoresis

Conforms

Single band of molecular weight > 200 KD

Reconstitution

Conforms

Each vial when reconstituted with 100 microliters of H₂O contains 2 micrograms of C. difficile toxin A enterotoxin in 50 mM Tris, pH 7.5, 50 mM Sodium Chloride, and 0.1% Trehalose.

Specification: PRD.0.ZQ5.10000038606

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.