


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
5-Formyluracil - 98%

Product Number: 414174
CAS Number: 1195-08-0
MDL: MFCD00192185
Formula: C₅H₄N₂O₃
Formula Weight: 140.1 g/mol



TEST

Specification

Appearance (Color)	White to Off-White
Appearance (Form)	Conforms to Requirements
Powder or Crystalline Powder	
Infrared spectrum	Conforms to Structure
Carbon	41.8 - 43.9 %
Nitrogen	19.5 - 20.5 %
Purity (HPLC)	≥ 97.5 %

Specification: PRD.0.ZQ5.10000011389

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.