


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Amino-2-methoxybenzoic acid - 96%

Product Number: 691429
CAS Number: 861306-04-9
MDL: MFCD09701221
Formula: C₈H₉NO₃
Formula Weight: 167.16 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)

Conforms to Requirements

White to Orange to Tan

Appearance (Form)

Powder or Crystals

Infrared spectrum

Conforms to Structure

Purity (Titration by NaOH)

95.5 - 104.5 %

Purity GC

≥ 95.5 %

Specification: PRD.0.ZQ5.10000015878

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.