


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sinapic acid - ≥98%, powder

Product Number: D7927
CAS Number: 530-59-6
MDL: MFCD00004401
Formula: C₁₁H₁₂O₅
Formula Weight: 224.21 g/mol



TEST

Specification

Appearance (Color)	Faint Yellow to Very Dark Yellow and Yellow-Brown
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Titration with NaOH	≥ 98 %
Purity (GC)	≥ 98 %
Recommended Retest Period	-----
5 years	

Specification: PRD.0.ZQ5.10000037753

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.