


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Methyl diphenylphosphinite - 97%

Product Number: 149497
CAS Number: 4020-99-9
MDL: MFCD00048025
Formula: C₁₃H₁₃OP
Formula Weight: 216.22 g/mol
Storage Temperature: 2 - 8 °C



| TEST | Specification |
|--------------------|-----------------------|
| Appearance (Color) | Colorless |
| Appearance (Form) | Liquid |
| Infrared spectrum | Conforms to Structure |
| Purity (GC) | ≥ 96.5 % |

Specification: PRD.0.ZQ5.10000040643

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.