

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

 β -Nicotinamide adenine dinucleotide, reduced disodium salt – ~98%, pkg of 0.2 mg (per vial)

Product Number:

N6660

CAS Number:

606-68-8

MDL:


MFCD00036200

Formula:

C₂₁H₂₇N₇Na₂O₁₄P₂

Formula Weight:

709.4 g/mol



TEST	Specification
------	---------------

Solubility (Color)

Colorless

0.2 mg plus 3 ml of 0.01 N NaOH

Solubility (Turbidity)

Clear

0.2 mg plus 3 ml of 0.01 N NaOH

Solubility (Form)

Solution

0.2 mg plus 3 ml of 0.01 N NaOH

mg/vial

0.2 - 0.3

Suitability

Suitable

Suitable for use as a co-factor in enzymatic assays. Gives comparable response as previous lots in an LDH assay system.

Specification: PRD.2.ZQ5.10000007205

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.