

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ambrette seed absolute - FGProduct Number: **W205001**
CAS Number: 8015-62-1
MDL: MFCD00239610

TEST	Specification
Appearance (Color)	Light Yellow to Amber
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.468 - 1.485
Infrared spectrum	Conforms to Structure
Specific Gravity	0.898 - 0.960
At 20°C	
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period	-----
5 Years	

Specification: PRD.2.ZQ5.1000008637

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.