


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Pyridinium dichromate - 98%

Product Number: 214698
CAS Number: 20039-37-6
MDL: MFCD00013105
Formula: C₁₀H₁₀N₂ · H₂Cr₂O₇
Formula Weight: 376.20 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Conforms to Requirements
Orange to Brown	
Appearance (Form)	Powder or Crystals or Powder with Chunk(s) or Crystals with Chunk(s)
Infrared spectrum	Conforms to Structure
Titration by Na ₂ S ₂ O ₃	26.9 - 28.3 %
% Chromium	

Specification: PRD.0.ZQ5.10000087990

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.