


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Thiopheneacetyl chloride – 98%

Product Number: **195995**
CAS Number: 39098-97-0
MDL: MFCD00005456
Formula: C₆H₅ClOS
Formula Weight: 160.62 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless to Brown
Appearance (Form)	Liquid
Proton NMR spectrum	Conforms to Structure
Purity GC VPCS	≥ 97.5 %

Specification: PRD.0.ZQ5.10000043677

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.