


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
tert-Butanol - TEBOL® 99, ≥99.3%

Product Number: **B85927**
CAS Number: 75-65-0
MDL: MFCD00004464
Formula: C₄H₁₀O
Formula Weight: 74.12 g/mol



TEST	Specification
Appearance (Color) White (Solid) or Colorless (Liquid)	Colorless or White
Appearance (Form) Solid or Liquid	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Supplier - Purity	≥ 99.30 %
Supplier - Water Content	≤ 0.15 %
Acidity as Acetic Acid Lyondell Data	≤ 0.0030
Color Test Lyondell Data	≤ 10 APHA
Non-volatile Matter G/100mL Lyondell Data	≤ 0.0010
Product of Supplier Lyondell	Confirmed

Specification: PRD.1.ZQ5.10000015891

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.