


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
3-Aminophenyl sulfone - 97%

**Product Number:** A74602  
**CAS Number:** 599-61-1  
**MDL:** MFCD00007792  
**Formula:** C<sub>12</sub>H<sub>12</sub>N<sub>2</sub>O<sub>2</sub>S  
**Formula Weight:** 248.3 g/mol



### TEST

### Specification

Appearance (Color) White to Off-White, Beige or Light Brown	Conforms to Requirements
Appearance (Form) Powder or Crystalline Powder and/or Chunks	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Carbon	56.0 - 60.1 %
Nitrogen	10.8 - 13.2 %
Sulfur (S)	12.4 - 13.4 %
Purity (HPLC)	≥ 96.5 %

Specification Date : 04/06/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.