

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Cellulose acetate propionate – average $M_n \sim 75,000$ by GPC, powder

Product Number:

330183

CAS Number:

9004-39-1

MDL:


MFCD00146341

Formula:

C₇₆H₁₁₄O₄₉

Formula Weight:

1,811.69 g/mol

**TEST****Specification**

Appearance (Color)

White

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Viscosity

15.0 - 25.0 s

Falling Ball, seconds

Water (by Karl Fischer)

≤ 3.00 %

Acetyl

≤ 2.50 %

Propionyl

45.00 - 49.00 %

Free Acid

≤ 0.03 %

Specification: PRD.0.ZQ5.10000051704

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.