

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Chorionic gonadotropin human - lyophilized powder, vial of ~10,000 IU

Product Number:

CG10

CAS Number:

9002-61-3

MDL:

MFCD00130797

Storage Temperature:

-20 °C

TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear to Slightly Hazy
1 vial + 10 mL, H ₂ O	
Potency	8000 - 12500
IU/Vial (based on the USP method for human chorionic gonadotropin)	
Tested For Infectious Agents	Tested

Specification: PRD.3.ZQ5.10000013308

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.