

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Platinum on carbon - Evonik Noblyst® P8075 1% Pt

Product Number:

907820

Pt

TEST	Specification
Appearance (Color)	Black
Appearance (Form)	Powder
ICP Major Analysis Confirms Platinum Component	Confirmed
Loss on Drying Calculated to dry basis	45.0 - 65.0 % wt
Platinum (Pt) Calculated to dry basis	0.9 - 1.1 % wt
Activity Crotonic acid hydrogenation ml/min.g	≥ 60
Product of Supplier Evonik	Confirmed
Registered Trademark Noblyst is a registered trademark of Evonik Degussa GmbH	Note
Recommended Retest Period 2 Years	-----

Specification: PRD.0.ZQ5.10000100619

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.