


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Potassium 5-chloropyridine-3-trifluoroborate - 97%

Product Number: **728845**



Formula: C₅H₃BClF₃KN
Formula Weight: 219.44 g/mol

TEST	Specification
Appearance (Color)	White and Faint to Dark Yellow and Faint Brown to Dark Brown
Appearance (Form)	Powder or Crystals
Proton NMR spectrum	Conforms to Structure
Fluorine NMR Spectra	Conforms to Structure
Carbon	26.4 - 28.3 %
Nitrogen	6.2 - 6.6 %
Purity (HPLC)	≥ 96.5 %

Specification: PRD.1.ZQ5.10000007389

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.