

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
ANTI-FLAG® M1 Agarose Affinity Gel

Product Number: **A4596**

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Ligand Immobilized	≥ 2.7
mg protein/ml packed gel	
Appearance (Color)	White
Appearance (Form)	Suspension
Gel Content	1.7 - 2.1 ml
1.7 - 2.1 ml resin slurry will give 1 ml PGV	
Binding Capacity / ml PGV	≥ 0.6 mg
Purity	≥ 95 %

Specification: PRD.0.ZQ5.10000067799

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.