

3050 Spruce Street, Saint Louis, MO 63103, USA  
Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)  
Email USA: [techserv@sial.com](mailto:techserv@sial.com)  
Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

Product Specification

Product Name:  
Anti-Rabbit IgG (whole molecule) antibody produced in goat – whole antiserum, lyophilized powder

Product Number: R0881

TEST	Specification
Appearance (Form)	Powder
Appearance (Color)	Faint Beige to Light Beige
Solubility (Solvent)	0.01M Phosphate Buffered Saline pH 7.4
Solubility (Volume)	5 - 5 ml
for 250 tests/vial	
Solubility (Volume)	20 - 20 ml
for 1000 tests/vial	
Solubility (Volume)	100 - 100 ml
for 5000 tests/vial	
Solubility (Color)	Very Light Brown to Very Dark Brown and Dark Brown-Red to Very Dark Brown-Red
Solubility (Turbidity)	Clear to Slightly Hazy
Antibody Concentration mg/test	0.22 - 0.26
by percipitin test (PPT)	
Suitable for RIA	Yes
Recommended Retest Period	-----
5 Years	

Specification: PRE.1.ZQ5.10000015317

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.