

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-Bovine IgG (whole molecule) antibody produced in rabbit - whole antiserum

Product Number:

**B8395**

Storage Temperature:

-20 °C

**TEST****Specification**

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Very Light Brown to Very Dark Brown and Dark Brown-Red to Very Dark Brown-Red
Protein by Biuret	40 - 100 mg/ml
Antibody Concentration(Method)	Quantitative Precipitin Assay
Antibody Concentration	4.0 - 7.0 mg/ml
ELISA - Titer	1:150,000-1:300,000
Recommended Retest Period	-----
5 Years	

Specification: PRE.0.ZQ5.10000020904

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.