

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Oleoyl-L- α -lysophosphatidic acid sodium salt - $\geq 98\%$, solid

Product Number: **L7260**
 CAS Number: 22556-62-3
 MDL: MFCD00133427
 Formula: $C_{21}H_{41}O_7P \cdot xNa^+$
 Formula Weight: 436.52 g/mol
 Storage Temperature: -20 °C



TEST

Specification

Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Very Faint Yellow
Solubility (Turbidity)	Clear
10 mg/ml in CHCl ₃ :MeOH:Acetic Acid (95:5:5)	
Phosphorus (P)	6.2 - 7.3 %
Free acid basis	
Sodium (Na)	3.8 - 10.6 %
Purity (TLC)	≥ 98 %
Purity (GC)	≥ 98 %

Note

Formula and formula weight are for the free acid

Specification: PRD.1.ZQ5.10000007013

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.