SIGMA-ALDRICH®

sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103, USA Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name: Luria broth base, Miller's modified - Tryptone, yeast extract and sodium chloride

Product Number:

L1900

TEST	Specification
Appearance (Color)	Faint Yellow to Light Yellow to Light Brown
Appearance (Form)	Powder
Solubility (Color)	Faint Brown-Yellow to Brown-Yellow
Solubility (Turbidity)	Clear to Slightly Opalescent
1.55%, H2O	
pН	6.8 - 7.4
(1.55% in Water at Room Temperature)	
Microbiological Test	Pass
Expiration Dated Product	

Specification: PRD.5.ZQ5.1000008705

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.