

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Fibrin from human plasma – insoluble powder

Product Number:

F5386

CAS Number:

9001-31-4

MDL:

MFCD00131058

TEST	Specification
------	---------------

Appearance (Color)

White to Yellow

Appearance (Form)

Powder

Solubility (Color)

Faint Yellow to Dark Yellow

Solubility (Turbidity)

Turbid

25 mg + 1 mL 1 N NaOH

Nitrogen

12 - 16 %

Tested For Infectious Agents

Tested

Specification: PRD.1.ZQ5.10000013051

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.