


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,2,3,4-Tetrahydroquinoxaline - 95%

Product Number: 702366
CAS Number: 3476-89-9
Formula: C₈H₁₀N₂
Formula Weight: 134.18 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Solid, Powder, Crystals and/or Chunks
Infrared spectrum
Purity (GC)

White to Yellow and White to Orange
Conforms to Requirements
Conforms to Structure
≥ 94.5 %

Specification: PRD.2.ZQ5.10000038667

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.