


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Chymostatin - microbial

Product Number: C7268
CAS Number: 9076-44-2
MDL: MFCD00071059



Storage Temperature: -20 °C

TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Yellow
Solubility (Turbidity)	Clear
10 mg/mL, Glacial HOAc	
% Peptide content	≥ 70

Specification: PRD.2.ZQ5.10000005141

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.