


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(S,R,S)-AHPC-PEG5-NH₂ hydrochloride

Product Number: 901850



Formula: C₃₄H₅₄ClN₅O₉S
Formula Weight: 744.34 g/mol
Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Color)	White to Light Yellow
Appearance (Form)	Powder or Crystals or Solid or Chunks
¹ H NMR Spectrum	Conforms to Structure
Purity (HPLC)	≥ 90.0 %

Specification: PRD.0.ZQ5.10000101482

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.